



Contribution ID: 228

Type: Young Crystallographers Lightning Talks(+poster)

## XRPD and TEM as tools to determine crystallite size of nanocrystalline iron

*Wednesday, 16 March 2022 18:25 (20 minutes)*

Analysis of the crystallite size data with different approaches utilising X-ray diffraction might lead to ambiguous results. For this reason, the use of different analytical techniques is required to confirm the XRD analysis. In the presented study a nanocrystalline iron was also examined with transmission electron microscopy to compare complementary measurements of crystallite sizes.

**Primary author:** ALBRECHT, Aleksander (West Pomeranian University of Technology in Szczecin)

**Co-author:** Mr MARCIN, Sadłowski (In the presented study a nanocrystalline iron was examined with different approaches utilising X-ray diffraction and transmission electron microscopy.)

**Presenter:** ALBRECHT, Aleksander (West Pomeranian University of Technology in Szczecin)

**Session Classification:** Postersession

**Track Classification:** Young crystallographers Lightning Talks