

MLZ internal Sample Environment Meeting

News from TUM Sample Environment

Andreas Buchner, Peter Biber, Manuel Suarez





Index

News from TUM Sample Environment


- LT news
- HT news
- Sample Changer for LT systems news
- Participation in developments and experiments
- Manpower. The status up to now.

Low Temperature News

- **Adiabatic Demagnetization Refrigerator (ADR)**



Temperature



- + Reception and commissioning
- + Cryostat performance characterization (in progress)



TBD:

- Sample holders, integration with the control system, auxiliary services and mechanical interface



Low Temperature News


- 1 K Cryostat



High Temperature News

- Advanced Materials Laserofen-Steuerung

AdvMat Steuerung Laserofen




Modular design, suitable for Laser + HTF heating methods.

Material Science Group, Sample Environment and Hereon

Testing machine for compression, tensile and fatigue tests




Sample changers for LT



• Automatic Linear sample exchange @LT

- Linear sample changer for Standard CCR



- Automated sample changer for low temperatures for the standard 80 mm CCR
- Allows the measurement of up to 10 samples without operator attendance.
- **Instrument control group is working on the control system**

• Robotic arm for sample changer @LT

- Automatic sample changer for Compact Cryostat and Standard CCR

On this day we must give the discussion about SE + Robotics

From SE we have some comments or points to contribute to the discussion:

- At SE we do not plan to develop “robotics arms” for sample exchangers in LT, but,
- We will be active collaborators and we´ll support its implementation
- Inst. Control group / Robotics group (what is robot?) / Sample Environment/Instruments
- ...


Participation in developments and experiments

- **New Pressure Cell (SANS-1)**

- Higher operational P-range, up 2kbar
- Change in the cooling circuit, thereby reducing the gradient on the sample from 1°C (old version) to 0,05°C
- Larger opening angle
- Additional optical window


A: Thermisch-stationäre Analyse
 Temperatur
 Typ: Temperatur
 Einheit: °C
 Zeit: 1 s
 21.07.2023 15:51

60,288 Max
 60
 59,5
 59
 58,5
 58
 57,5
 56,2
 39,1
 22 Min




A: Thermisch-stationäre Analyse
 Temperatur
 Typ: Temperatur
 Einheit: °C
 Zeit: 1 s
 21.07.2023 17:03

60,022 Max
 58,2
 58,1
 58
 57,9
 50,72
 43,54
 36,36
 29,18
 22 Min



Manpower

- Manuel Suarez (Management and LT)
- Andreas Buchner (HP, HT)
- Peter Biber (LT)
- Jank Wenzlaff (Workshop)
- Edith Rutenfranz (Student)



Stellenangebote Forschungs-Neutronenquelle
Heinz Maier-Leibnitz (FRM II)

+ 2 positions for Engineer
+ 1 position for Scientist



Thanks you for your time

Manuel Suarez

Manuel.Suarez@frm2.tum.de

+49 (0)89 289-10769

Andreas Buchner, Peter Biber, Jank Wenzlaff and Edith Rutenfranz

Thank you for your time



MLZ is a cooperation between:

Sample Environment @ TUM Summary

High Temperature ~ (6 HTF, 4 IRF):
RT-1900K

Low Temperature ~ (19CCR, 6He3, 2DU):
RT- 50mK

Magnetic Field: (~5MG)
Up 12T

High Pressure: ~(5HPC)
Up 7,5kbar

